Name Date

Period

Grade:

EXPERIMENT 38 ACID-BASE TITRATION

PRELAB QUESTIONS:

- 1. Define these terms: standard solution, titration, endpoint, neutralization reaction.
- 2. What color does the indicator phenolphthalein turn in base solutions? What color does it turn in acid solutions?
- 3. Write a balance equation for the following neutralization reactions:
 - a) $HBr + KOH \rightarrow$
 - b) $HI + Ca(OH) \rightarrow$
 - c) $H_2SO_4 + NaOH \rightarrow$
 - d) $H_2SO_4 + Mg(OH)_2 \rightarrow$
- 4. If 30.0 mL of 0.5 M KOH are needed to neutralize 10.0 mL of HBr of unknown concentration, what is the molarity of the HBr?
- 5. What volume of 0.3 M HI is needed to exactly neutralize 50.0 mL of 0.3 M KOH?
- 6. What is the molarity of the base if 65.0 mL of 0.1 M H₂SO₄ is neutralized by 35.0 mL of Ba(OH)₂ solution ?
- 7. It took 150.0 mL Mg(OH)₂ to neutralize 75.0 mL of H₂SO₃. What is the molarity of the H₂SO₃?

Name Period

Date Lab Partners

EXPERIMENT 38 ACID-BASE TITRATION

OBSERVATIONS AND DATA:

DATA TABLE

| DITTI TIBEE | | | | | | | | | | | |
|-----------------|---------|------|---------|------|---------|------|---------|------|--|--|--|
| | Trial 1 | | Trial 2 | | Trial 3 | | Trial 4 | | | | |
| | HCI | NaOH | HCI | NaOH | HCI | NaOH | HCI | NaOH | | | |
| Initial reading | | | | | | | | | | | |
| Final reading | | | | | | | | | | | |
| Volume used | | | | | | | | | | | |

CALCULATIONS:

For each trial calculate the molarity of the NaOH solution using the relationship:

 $\mathbf{M}_{\mathbf{a}} \times \mathbf{V}_{\mathbf{a}} = \mathbf{M}_{\mathbf{b}} \times \mathbf{V}_{\mathbf{b}}$

Show all work!

Trial 1

Trial 2

Trial 3

Trial 4

CONCLUSIONS AND QUESTIONS:

- 1. List possible sources of error in this experiment.
- 2. How reproducible were the results of the four trials? How did your results compare with those of the other lab groups?
- 3. How many mL of 0.200 M LiOH are needed to titrate 50.0 mL of 0.500 M H₂SO₄?

Discussion

Conclusion